

Genotyping Protocol: **MMRRC 203**

**Assay Type:** PCR- can distinguish heterozygous animals from homozygous animals

**DNA Extraction:** DNA from tail snips was extracted using Qiagen's DNeasy Blood and Tissue kit (Cat# 69506). Kit directions for animal tissues were performed with a few minor modifications as follows: repeat AW1 and AW2 wash steps one time, elute in 200µl of AE buffer once.

**Strain Description:** The 3<sup>rd</sup> and 4<sup>th</sup> exons of *Sod1* were replaced with a promoter-less neomycin resistance gene (CITEneo). Details can be found in Huang, TT *et al.* (1997) Arch Biochem Biophys 344(2):424-32.




Image from Huang, TT *et al.* (1997) Arch Biochem Biophys 344(2):424-32.

**Primer Information:**

- |                     |   |
|---------------------|---|
| 1) Name: M203 Mut F | Sequence: 5'-CCG TCT TTT GGC AAT GTG AGG-3' |
| 2) Name: M203 Mut R | Sequence: 5'-TGG GGT ACC TTC TGG GCA TCC-3' |
| 3) Name: M203 WT R  | Sequence: 5'-GGC GGA TGA AGA GAG GTG AGC-3' |
| 4) Name: M203 WT R  | Sequence: 5'-ATT GGC CAC ACC GTC CTT TCC-3' |

**Primer location:** M203 Mut F and M203 Mut R are located in the CITEneo construct. M203 WT F is located in Exon 3 of *Sod1* and M203 WT R is located in Exon 4 of *Sod1*.

**Assay name: Sod1 KO PCR****PCR Master Mix Components:**

component	manufacturer	concentration	µl/rxn
Buffer with MgCl <sub>2</sub> (green cap)	Roche	10X	2
dNTPs	Promega (Cat# U1515)	1.25mM	3.2
M203 Mut F	Sigma or IDT	25µM	0.3
M203 Mut R	Sigma or IDT	25µM	0.3
M203 WT F	Sigma or IDT	25µM	0.3
M203 WT R	Sigma or IDT	25µM	0.3
FastStart Taq	Roche (Cat# 12032953001)	5 U/µl	0.2
sterile water			12.4

**PCR Setup:**

Final Reaction: 19µl master mix & 1µl DNA template (10-20ng/µl)

All reactions were performed in 200µl thin walled PCR tubes and were run in Perkin Elmer 2400 thermocycler or Applied Biosystems 2700 thermocycler.

**Cycle Parameters:**

- 1) 95°C 4 minutes
- 2) 94°C 30 seconds
- 3) 66°C 30 seconds
- 4) 72°C 1 minute
- 5) Repeat steps 2-4 34 times for a total of 35 cycles
- 6) 72°C 7 minutes
- 7) 4°C hold until refrigerate product

**Product Analysis:**

All products were analyzed on the Qiaxcel (instrument and all supplies from Qiagen) with the Qiaxcel DNA Screening Kit (Cat# 929004).

Alignment Marker: QX Alignment Marker 15bp/1kb (Cat# 929521)

Size Marker: QX DNA Size Marker 50bp-800bp (Cat# 929556)

Method: AM320

Injection: 10s at 5KV


Separation: 320s at 6KV

**Expected products:**

Wild type allele: 780bp product

Mutant allele: 350bp product

**Example gel:**



Lane A12 displays a wild-type sample (780bp product)

Lane B1 displays a heterozygous sample (780bp and 350bp products)

\*Please note: the 15bp and 3 Kb bands are reference markers for the QIAxcel method and do not represent expected products.\*