

Genotyping Protocol: **MMRRC line 30088**

**Assay Type:** PCR - can distinguish heterozygous animals from homozygous animals

**DNA Extraction:** DNA from tail snips was extracted using Sigma's RedExtract-N-Amp Tissue PCR Kit (Cat#XNAT2R). Kit directions for animal tissues were performed with a few minor modifications as follows: Use only 50 µl of Extraction Solution, 12.5 µl Tissue Preparation Solution and 50 µl of Neutralization Solution B.

**Mutation Information:** In this line most of exon 1 of the *Zip4* gene (*Slc39a4*, a gene involved in zinc absorption) has been replaced with EGFP. EGFP was inserted in-frame.

**Primer Information:**

- |                      |   |
|----------------------|---|
| 1) Name: mZip4WT(as) | Sequence: 5'-GAA TGG AGC TGT GTG TCC CAG G-3'       |
| 2) Name: mZip4WT(s)  | Sequence: 5'-CAC AGG GAC TTG TGT TGG CTG TG-3'      |
| 3) Name: mZip4Mut(s) | Sequence: 5'-TCA CTG CAT TCT AGT TGT GGT TTG TCC-3' |

**Primer location:** mZip4WT(as) binds to intron 1 of *Zip4*. mZip4WT(s) binds to exon 1 of *Zip4*. mZip4Mut(s) binds to the EGFP open reading frame.

**Assay Name:** Zip4 KO PCR

**PCR Master Mix Components:**

component	manufacturer	concentration	µl/rxn
Extract-N-Amp PCR Reaction Mix	Sigma (Cat#XNAT2R)	2X	10
mZip4WT(as)	IDT or Sigma	25µM	0.3
mZip4WT(s)	IDT or Sigma	25µM	0.3
mZip4Mut(s)	IDT or Sigma	25µM	0.3
sterile water			5.1

**PCR Setup:**

Final Reaction: 16µl master mix & 4µl DNA template (10-20 ng/µl)

All reactions were performed in 200µl thin walled PCR tubes and were run in Perkin Elmer 2400 thermocycler or Applied Biosystems 2700 thermocycler.

**Cycle Parameters:**

- 1) 94°C      3 minutes
- 2) 94°C      1 minute
- 3) 66°C      1 minute
- 4) 72°C      1 minute
- 5) Repeat steps 2-4    34 times for a total of 35 cycles
- 6) 72°C      10 minutes
- 7) 4°C      hold until refrigerate product

**Product Analysis:**


All products were analyzed on a 3% agarose gel with ethidium bromide staining

**Expected products:**

Wild-type: 375 bp  
Heterozygous: 375 bp and 277 bp  
Homozygous: 277 bp

Control DNA: positive and negative animals

**Example of Gel:**



Wells 1-4 are heterozygous. Well 7 is a wild-type control. Well 8 is 1 Kb Plus DNA ladder (Invitrogen Cat. # 10787-018).