


Genotyping Protocol: **MMRRC 30397**

Assay Type: PCR- can distinguish heterozygous animals from homozygous animals

DNA Extraction: DNA from tail snips was extracted using Qiagen's DNeasy Blood and Tissue kit (Cat# 69506). Kit directions for animal tissues were performed with a few minor modifications as follows: repeat AW1 and AW2 wash steps one time, elute in 200µl of AE buffer once.

Strain Description: This strain has IRES-GFP^{LacZ} inserted into exon 1 of the reticulon 4 (*Rtn4*) gene on Chromosome 11, resulting in the deletion of exon 1 downstream of the ATG start codon. This abolishes *Nogo-A* and *Nogo-B* expression. Details can be found in Zheng et al (2003) Neuron 38:213-24.



Primer Information:

- | | |
|--------------------------|--|
| 1) Name: M30399 NNGenF2 | Sequence: 5'-CAG TAG CTG CAG CAT CAT CG-3' |
| 2) Name: M30399 NNGenR3 | Sequence: 5'-CTC TCC AGC ACC TCC AAT TC-3' |
| 3) Name: M30399 IRES200R | Sequence: 5'-AGA GGA ACT GCT TCC TTC AC-3' |

Primer location: M30399 NNGenF2 and NNGenR3 are located in exon 1 of the *Rtn4* gene; M30399 IRES200R is located in the inserted IRES-GFP^{LacZ}.

Assay name: Nogo Deletion (B1) PCR

Mutant PCR:

PCR Master Mix Components:

component	manufacturer	concentration	µl/rxn
Buffer with MgCl ₂ (green cap)	Roche	10X	2
dNTPs	Promega (Cat# U1515)	1.25mM	3.2
M30399 NNGenF2	Sigma	25µM	0.3
M30399 IRES200R	Sigma	25µM	0.3
FastStart Taq	Roche (Cat# 12032953001)	5 U/µl	0.2
sterile water			13

PCR Setup:

Final Reaction: 19µl master mix & 1µl DNA template (10-20ng/µl)

10.06.08 MS
07.30.10 HB updated

All reactions were performed in 200µl thin walled PCR tubes and were run in Perkin Elmer 2400 thermocycler or Applied Biosystems 2700 thermocycler.

Cycle Parameters:


- 1) 95°C 5 minutes
- 2) 94°C 1 minute
- 3) 68°C 1 minute
- 4) 72°C 1 minute
- 5) Repeat steps 2-4 34 times for a total of 35 cycles
- 6) 72°C 10 minutes
- 7) 4°C hold until refrigerate product

Product Analysis:

All products were analyzed on a 3% agarose gel with ethidium bromide staining.

Expected products: 350bp mutant band

Example gel:



Wells 1 and 2 are positive for the mutant allele. Well 3 and 4 are blanks, and Well 5 is a wild-type control, and well 6 is a heterozygous control. Well 7 is 1Kb+ Ladder (Invitrogen Cat# 10787-018).

WT PCR:

PCR Master Mix Components:

component	manufacturer	concentration	µl/rxn
Buffer with MgCl ₂ (green cap)	Roche	10X	2
dNTPs	Promega (Cat# U1515)	1.25mM	3.2
M30399 NNGenF2	Sigma	25µM	0.3
M30399 NNGenR3	Sigma	25µM	0.3
FastStart Tag	Roche (Cat# 12032953001)	5 U/µl	0.2
sterile water			13

PCR Setup:

Final Reaction: 19µl master mix & 1µl DNA template (10-20ng/µl)

All reactions were performed in 200µl thin walled PCR tubes and were run in Perkin Elmer 2400 thermocycler or Applied Biosystems 2700 thermocycler.

10.06.08 MS
07.30.10 HB updated

Cycle Parameters:


- | | | |
|----|--|--------------------------------|
| 1) | 95°C | 5 minutes |
| 2) | 94°C | 1 minute |
| 3) | 62°C | 1 minute |
| 4) | 72°C | 1 minute |
| 5) | Repeat steps 2-4 34 times for a total of 35 cycles | |
| 6) | 72°C | 10 minutes |
| 7) | 4°C | hold until refrigerate product |

Product Analysis:

All products were analyzed on a 3% agarose gel with ethidium bromide staining.

Expected products: 200bp wild type band

Example gel:



Wells 10 and 11 are positive for the WT allele. Wells 12 and 13 are blanks; Well 14 is a WT control, and Well 15 is a heterozygous control. Well 16 is 1Kb+ Ladder (Invitrogen Cat# 10787-018).

Product Analysis:

Wild Type: 200bp on WT gel, no product on Mut gel

Heterozygous: 200bp on WT gel, 350bp on Mut gel

Homozygous: no product on WT gel, 350bp on Mut gel