


Genotyping Protocol: **MMRRC 30665**

Assay Type: PCR- can distinguish heterozygous animals from homozygous animals

DNA Extraction: DNA from tail snips was extracted using Qiagen's DNeasy Blood and Tissue kit (Cat# 69506). Kit directions for animal tissues were performed with a few minor modifications as follows: repeat AW1 and AW2 wash steps one time, elute in 200ul of AE buffer once.

Strain Characteristics: This strain has the LacZ gene knocked in to the T-cell acute lymphocytic leukemia 1 gene (*Tal1*), which results in the deletion of exons 4 and 5. This strain is homozygous lethal at ~day 8.5 due to failure of vascular system development in the embryo. Details can be found in Elefanty et al (1998) Proc Natl Acad Sci USA 95(20):11897-902.



Primer Information:

- | | |
|--------------------|---|
| 1) Name: Tal1.mutF | Sequence: 5'-ATA TTG CTG AAG AGC TTG GCG GC-3' |
| 2) Name: Tal1.wtF | Sequence: 5'-GTT TTG GTC TAG AGT TTG TGA GCC-3' |
| 3) Name: Tal1.R | Sequence: 5'-GCA TGC TCA AGG CTG CTG ACT TGG-3' |

Primer location: Tal1.mutF is located in the inserted LacZ cassette. Tal1.wtF and Tal1.R are located in exon 5.

Assay name: Tal1 KI PCR

PCR Master Mix Components:

component	manufacturer	concentration	µl/rxn
Buffer with MgCl ₂ (green cap)	Roche	10X	2
dNTPs	Promega (Cat# U1515)	1.25mM	3.2
Tal1.mutF	Sigma	25µM	0.3
Tal1.wtF	Sigma	25µM	0.3
Tal1.R	Sigma	25µM	0.3
FastStart Taq	Roche (Cat# 12032953001)	5 U/µl	0.2
sterile water			12.7

PCR Setup:

Final Reaction: 19µl master mix & 1µl DNA template (10-20ng/µl)

All reactions were performed in 200µl thin walled PCR tubes and were run in Perkin Elmer 2400 thermocycler or Applied Biosystems 2700 thermocycler.

03.15.11 MS
02.24.14 MLS

Cycle Parameters:

- 1) 95°C 3 minutes
- 2) 94°C 30 seconds
- 3) 66°C 30 seconds
- 4) 72°C 30 seconds
- 5) Repeat steps 2-4 34 times for a total of 35 cycles
- 6) 72°C 10 minutes
- 7) 4°C hold until refrigerate product

Product Analysis:

All products were analyzed on the Qiaxcel (instrument and all supplies from Qiagen) with the Qiaxcel DNA Screening Kit (Cat# 929004).

Alignment Marker: QX Alignment Marker 15bp/1kb (Cat# 929521)

Size Marker: QX DNA Size Marker 50bp-800bp (Cat# 929556)

Method: AM320 Injection: 10s at 5KV


Separation: 320s at 6KV

Expected products:

WT: 424bp

Mut: ~650bp

Example gel:



Lane B05 displays a WT sample (424bp band).

Lane B06 displays a heterozygous sample (424bp and 650bp bands).

Please note: the 15bp and 1kb bands are reference markers specific to the QIAxcel method and do not represent expected products.