


Genotyping Protocol: **MMRRC 30748**

Assay Type: PCR – can distinguish heterozygous animals from homozygous animals

DNA Extraction: DNA from tail snips was extracted using Qiagen's DNeasy Blood and Tissue kit (Cat# 69506). Kit directions for animal tissues were performed with a few minor modifications as follows: repeat AW1 and AW2 wash steps one time, elute in 200µl of AE buffer once.

Strain Description: This strain carries an insertion of a PLAP secretory trap cassette within the semaphorin 4b gene (*Sema4b*) on Chromosome 7. The PLAP secretory trap vector contains β -geo, PLAP(human placental alkaline phosphatase), and an internal ribosome entry site (IRES). Details can be found in Friedel et al (2005) Proc Natl Acad Sci U S A 102:13188-93.



Primer Information:

- | | |
|--------------------------|--------------------------------------|
| 1) Name: M30748 Mut IRES | Sequence: TGC ACA TGC TTT ACG TGT G |
| 2) Name: M30748 Mut PLAP | Sequence: TGC CGC GTG TCG TGT TGC AC |
| 3) Name: M30748 WT F | Sequence: AAT CTG GAT CAG TGG ACA GC |
| 4) Name: M30748 WT R | Sequence: CCA CAA GCT GCT GGA ATC CT |

Primer location: M30748 WT F is located in intron 14 and M30748 WT R is located in exon 15 of *Sema4b*. M30748 Mut IRES and Mut PLAP are located in the inserted cassette.

Assay Name: *Sema4b* PCR

Mut PCR:

PCR Master Mix Components:

Component	manufacturer	concentration	µl/rxn
Buffer with MgCl ₂ (green cap)	Roche	10X	2
dNTP	Promega (Cat# U1515)	1.25mM	3.2
M30748 Mut IRES	Sigma	25µM	0.3
M30748 Mut PLAP	Sigma	25µM	0.3

11.16.09 MS
08.02.10 HB updated

FastStart Taq	Roche (Cat# 12032953001)	5 U/μl	0.2
sterile water			13

PCR Setup:

Final Reaction: 19μl master mix & 1μl DNA template (10-20ng/μl)

All reactions were performed in 200μl thin walled PCR tubes and were run in Perkin Elmer 2400 thermocycler or Applied Biosystems 2700 thermocycler.

Cycle Parameters:

- 1) 95°C 3 minutes
- 2) 94°C 30 seconds
- 3) 55°C 30 seconds
- 4) 72°C 30 seconds
- 5) Repeat steps 2-3 34 times for a total of 35 cycles
- 6) 72°C 10minutes
- 7) 4°C hold until refrigerate product

Product Analysis:

All products were analyzed on the Qiaxcel (instrument and all supplies from Qiagen) with the Qiaxcel DNA Screening Kit (Cat# 929004).

Alignment Marker: QX Alignment Marker 15bp/1000bp (Cat# 929521)

Size Marker: QX DNA Size Marker 50-800bp (Cat# 929556)

Method: AH320 Injection: 20s at 2kV

Separation: 320s at 6kV

Expected Product: 510bp

WT PCR:

PCR Master Mix Components:

Component	manufacturer	concentration	μl/rxn
Buffer	Roche	10X	2
dNTP	Promega (Cat# U1515)	1.25mM	3.2
M30748 WT F	Sigma	25μM	0.3
M30748 WT R	Sigma	25μM	0.3
FastStart Taq	Roche (Cat# 12032953001)	5 U/μl	0.2
sterile water			13

PCR Setup:

Final Reaction: 19μl master mix & 1μl DNA template (10-20ng/μl)

All reactions were performed in 200μl thin walled PCR tubes and were run in Perkin Elmer 2400 thermocycler or Applied Biosystems 2700 thermocycler.

Cycle Parameters:

- 1) 95°C 3 minutes
- 2) 94°C 30 seconds
- 3) 66°C 30 seconds
- 4) 72°C 30 seconds
- 5) Repeat steps 2-3 34 times for a total of 35 cycles
- 6) 72°C 10minutes
- 7) 4°C hold until refrigerate product

11.16.09 MS
08.02.10 HB updated

Product Analysis:

All products were analyzed on the Qiaxcel (instrument and all supplies from Qiagen) with the Qiaxcel DNA Screening Kit (Cat# 929004).

Alignment Marker: QX Alignment Marker 15bp/1000bp (Cat# 929521)

Size Marker: QX DNA Size Marker 50-800bp (Cat# 929556)

Method: AH320


Injection: 20s at 2kV

Separation: 320s at 6kV

Expected Product: 490bp

Genotype	Mut PCR	WT PCR
WT	no product	490bp
Heterozygous	510bp	490bp
Homozygous	510bp	no product

Example Gel:



Lane C02 displays a sample positive for the mutant allele (510bp band).

Lane C08 displays a sample positive for the WT allele (490bp band).

Please note: the bands at 15bp and 1000bp are reference bands specific to the QIAxcel method and do not represent expected products.