

11.11.09 MS  
 08.02.10 HB updated  
 03.04.14 MLS

## Genotyping Protocol: **MMRRC 30963**

### **Current Background Strain of M30963: C57BL/6**


Details on the importance of background strain can be found in Xiao et al (2009) *Biol Blood Marrow Transplant* 15:1-11.

*\*This strain has an identical genetic alteration to M30964 – the two strains differ in current background strain.\**

**Assay Type:** PCR- can distinguish heterozygous animals from homozygous animals

**DNA Extraction:** DNA from tail snips was extracted using Qiagen's DNeasy Blood and Tissue kit (Cat# 69506). Kit directions for animal tissues were performed with a few minor modifications as follows: repeat AW1 and AW2 wash steps one time, elute in 200ul of AE buffer once.

**Strain Description:** This strain has a targeting construct inserted into the DNA cross-link repair 1C, PSO2 homolog gene (*Dclre1c*) on Chromosome 2, resulting in the deletion of exons 6-11 of *Dclre1c*. The targeting construct is made up of a 1.8kb Neomycin gene flanked by exons 4-6 and 11-14 of *Dclre1c*. Details can be found in Li et al (2005) *J. Immunol.* 174:2420-2428.



**Primer Information:**

- |                       |   |
|-----------------------|---|
| 1) Name: M30964 E6 F  | Sequence: 5'-TCA GGG CAG TAA TGG AAC TGT C-3' |
| 2) Name: M30964 E6 R  | Sequence: 5'-CTG TGT GCT CAC ACA TGC AC-3'    |
| 3) Name: M30964 Neo R | Sequence: 5'-TGC CTG CTC TTT ACT GAA GG-3'    |

**Primer location:** M30964 E6 F is located in exon 6, and M30964 E6 R is located in intron 6 of the *Dclre1c* gene on Chromosome 2. M30964 Neo R is located in the targeting construct.

**Assay name: Dclre1c PCR**

**PCR Master Mix Components:**

component	manufacturer	concentration	μl/rxn
Buffer with MgCl <sub>2</sub> (green cap)	Roche	10X	2
dNTP	Promega (Cat# U1515)	1.25mM	3.2
M30964 E6 F	Sigma	25μM	0.3
M30964 E6 R	Sigma	25μM	0.3
M30964 Neo R	Sigma	25μM	0.3
FastStart <i>Taq</i>	Roche (Cat# 12032953001)	5 U/μl	0.2
sterile water			12.7

11.11.09 MS  
08.02.10 HB updated  
03.04.14 MLS

**PCR Setup:**

Final Reaction: 19µl master mix & 1µl DNA template (10-20ng/µl)

All reactions were performed in 200µl thin walled PCR tubes and were run in Perkin Elmer 2400 thermocycler or Applied Biosystems 2700 thermocycler.

**Cycle Parameters:**

- 1) 95°C 5 minutes
- 2) 94°C 30 seconds
- 3) 63°C 30 seconds
- 4) 72°C 30 seconds
- 5) Repeat steps 2-4 34 times for a total of 35 cycles
- 6) 72°C 10 minutes
- 7) 4°C hold until refrigerate product

**Product Analysis:**

All products were analyzed on the Qiaxcel (instrument and all supplies from Qiagen) with the Qiaxcel DNA Screening Kit (Cat# 929004).

Alignment Marker: QX Alignment Marker 15bp/3Kb (Cat# 929522)

Size Marker: QX DNA Size Marker 100-3Kb (Cat# 929553)

Method: AM320

Injection: 10s at 5KV


Separation: 320s at 6KV

Heterozygous: 169bp, 400bp

Homozygous for mutant allele: 400bp

Wild Type: 169bp

**Example gel:**



Lane F07 displays a WT sample (169bp band).  
Lane F08 displays a heterozygous sample (169bp and 400bp bands).  
Lane F02 displays a homozygous mutant sample (400bp band).  
  
\*Please note: the 15bp and 3kb bands are reference markers specific to the QIAxcel method and do not represent expected products.\*