

05.16.11 MS  
01.07.14 MLS

## Genotyping Protocol: **MMRRC 32780**

**Assay Type:** PCR to detect transgene positive animals- cannot distinguish hemizygous animals from homozygous animals

**DNA Extraction:** DNA from tail snips was extracted using Qiagen's DNeasy Blood and Tissue kit (Cat# 69506). Kit directions for animal tissues were performed with a few minor modifications as follows: repeat AW1 and AW2 wash steps one time, elute in 200µl of AE buffer once.

**Strain Description:** This strain carries a transgene which contains Cre-ERT2 driven by the human doublecortin gene (*DCX*) promoter.

### Primer Information:

1) Name: hDCX.F	Sequence: 5'-TGA ATG TCG GAT AGC TGC AC-3'
2) Name: Cre.R	Sequence: 5'-GCA AAC GGA CAG AAG CAT TT-3'
3) Name: ERT2-F	Sequence: 5'-GAT TGG TCT CGT CTG GCG CTC C-3'
4) Name: ERT2-R	Sequence: 5'-ACG GCT AGT GGG CGC ATG T-3'

**Primer location:** hDCX.F is located in the promoter region of the human *DCX* gene which is contained in the transgene. Cre.R is located in the *Cre* gene. ERT2-F and ERT2-R bind to the human estrogen receptor T2 (ERT2) gene designed in the pCAG-ERT2-Cre-ERT2 vector.

**Assay name: M32780 PCR**

### hDCX-Cre PCR:

#### PCR Master Mix Components:

component	manufacturer	concentration	µl/rxn
Buffer with MgCl <sub>2</sub> (green cap)	Roche	10X	2
dNTPs	Promega (Cat# U1515)	1.25mM	3.2
hDCX.F	Sigma	25µM	0.3
Cre.R	Sigma	25µM	0.3
FastStart <i>Taq</i>	Roche (Cat# 12032953001)	5 U/µl	0.2
sterile water			13

### PCR Setup:

Final Reaction: 19µl master mix & 1µl DNA template (10-20ng/µl)

All reactions were performed in 200µl thin walled PCR tubes and were run in Perkin Elmer 2400 thermocycler or Applied Biosystems 2700 thermocycler.

### Cycle Parameters:

1)	95°C	3 minutes
2)	94°C	1 minute
3)	64°C	1 minute
4)	72°C	1 minute
5)	Repeat steps 2-4 34 times for a total of 35 cycles	
6)	72°C	10 minutes
7)	4°C	hold until refrigerate product

### Expected product:

Positive: 899bp

Negative: no product

05.16.11 MS  
01.07.14 MLS

**Product Analysis:**

All products were analyzed on the Qiaxcel (instrument and all supplies from Qiagen) with the Qiaxcel DNA Screening Kit (Cat# 929004).

Alignment Marker: QX Alignment Marker 15bp/3Kb (Cat# 929522)


Size Marker: QX DNA Size Marker 100-3Kb (Cat# 929553)

Method: AM320

Injection: 10s at 5KV

Separation: 320s at 6KV

**Example gel:**



Lane A09 displays a transgene negative sample (no product).  
Lane A10 displays a transgene positive sample (899bp product).

\*Please note: the bands appearing at 3kb and 15bp are reference bands specific to the Qiaxcel method and do not represent expected products.\*

**ERT2 PCR:**

**Note: it is not necessary to run this assay in addition to the first unless confirmation of the presence of the tamoxifen inducible element is desired.**

**PCR Master Mix Components:**

component	manufacturer	concentration	µl/rxn
Buffer with MgCl <sub>2</sub> (green cap)	Roche	10X	2
dNTPs	Promega (Cat# U1515)	1.25mM	3.2
ERT2-F	IDT	25µM	0.3
ERT2-R	IDT	25µM	0.3
FastStart Taq	Roche (Cat# 12032953001)	5 U/µl	0.2
sterile water			13

**PCR Setup:**

Final Reaction: 19 µl master mix & 1 µl extracted DNA (10-20ng/µl).

All reactions were performed in 200 µl thin walled PCR tubes and were run in Perkin Elmer 2400 thermocycler or Applied Biosystems 2700 thermocycler.

**Cycle Parameters:**

- 1) 95°C 3 minutes
- 2) 94°C 30 seconds
- 3) 60°C 30 seconds
- 4) 72°C 30 seconds
- 5) Repeat steps 2-4 34 times for a total of 35 cycles
- 6) 72°C 10 minutes
- 7) 4°C hold until refrigerate product

05.16.11 MS  
01.07.14 MLS

**Expected product:** 502 bp

**Product Analysis:**

Products may be analyzed on the Qiaxcel.

For analysis on the Qiaxcel (instrument and all supplies from Qiagen) with the Qiaxcel DNA Screening Kit (Cat# 929004):

Alignment Marker: QX Alignment Marker 15bp/3kb (Cat# 929522)


Size Marker: QX DNA Size Marker 100bp-3kb (Cat# 929553)

Method: AH320

Injection: 20s at 2kV

Separation: 320s at 6kV

**Example gel:**



<u>Lane</u>	<u>Sample</u>
A07	pCAG-ERT2-Cre-ERT2 (+)
B07	C57BL/6 (neg Control)

\*Please note: the bands appearing at 3kb and 15bp are reference bands specific to the Qiaxcel method and do not represent expected products.\*