

08.15.12 MLS


08.24.15 MLS

Genotyping Protocol: **MMRRC 36739**

Assay Type: PCR; Can distinguish heterozygous animals from homozygous animals.

DNA Extraction: DNA from tail snips was extracted using Qiagen's DNeasy Blood and Tissue kit (Cat# 69506). Kit directions for animal tissues were performed with a few minor modifications as follows: repeat AW1 and AW2 wash steps one time, elute in 200µl of AE buffer once.

Mutation Information: Exon 5 and part of intron 4 of the mouse catalase gene (*Cat*) was replaced by a neomycin cassette. Details can be found in Ho et al (2004) J Biol Chem 279(31):32804-12.



Primer Information:

- 1) Name: M36739 F Sequence: 5'-GGG ACG AAG ACT CCA GAA GTC CC-3'
- 2) Name: M36739 R Sequence: 5'-GCC TGG AGA ACA GGC TGT GCC-3'
- 2) Name: M36739 Neo Sequence: 5'-TTG GCG GCG AAT GGG CTG AC-3'

Primer Location: M36739 F is located in Exon 5 of *Cat*, and M36739 R is located in Intron 4 of the mouse *Cat* gene. M36739 Neo is located in the inserted neomycin cassette.

Assay Name: **Cat KO PCR**

PCR Master Mix Components:

component	manufacturer	concentration	µl/rxn
Buffer with MgCl ₂ (green cap)	Roche	10X	2
dNTPs	Promega (Cat# U1515)	1.25mM	3.2
M36739 F	Sigma	25µM	0.3
M36739 R	Sigma	25µM	0.3
M36739 Neo	Sigma	25µM	0.3
FastStart <i>Taq</i>	Roche (Cat# 12032953001)	5 U/µl	0.2
sterile water			12.7

PCR Setup:

Final Reaction: 19µl master mix & 1µl DNA template (10-20ng/µl)

All reactions were performed in 200µl thin walled PCR tubes and were run in an Applied Biosystems 2700 thermocycler.

08.15.12 MLS

08.24.15 MLS

Cycle Parameters:

- 1) 95°C 3 minutes
- 2) 94°C 30 seconds
- 3) 66°C 30 seconds
- 4) 72°C 1 minute
- 5) Repeat steps 2-4 34 times for a total of 35 cycles
- 6) 72°C 10 minutes
- 7) 4°C hold until refrigerate product

Product Analysis:

All products were analyzed on the Qiaxcel (instrument and all supplies from Qiagen) with the Qiaxcel DNA Screening Kit (Cat# 929004).

Alignment Marker: QX Alignment Marker 15bp/3Kb (Cat# 929522)

Size Marker: QX DNA Size Marker 100-3Kb (Cat# 929553)


Method: AM320

Injection: 10s at 5KV

Separation: 320s at 6KV

- Expected products:**
- Wild type: 450bp product only
 - Heterozygous: 450bp and 650bp products
 - Homozygous: 650bp product only

Example of Gel:



Lanes A1 and A2 display heterozygous samples (450bp and 650bp products).
Lanes A3 and A4 displays extraction and PCR blanks, respectively.
Lane A5 displays a wild type sample (450bp product only).
Lane A6 displays a homozygous mutant sample (650bp product only).

Please note: the 15bp and 3kb bands are reference markers specific to the Qiaxcel method and do not represent expected products.