# Genotyping Protocol: MMRRC 422

**Assay Type**: PCR - cannot distinguish hemizygous animals from homozygous animals, can only distinguish transgene positive from transgene negative animals.

**DNA Extraction**: DNA from tail snips was extracted using Sigma's Extract-N-Amp Tissue PCR Kit (Cat# XNAT2R). Kit directions for animal tissues were performed with a few minor modifications as follows: Use only 50 μl of Extraction Solution, 12.5 μl Tissue Preparation Solution and 50 μl of Neutralization Solution B.

#### **Primer Information:**

1) Name: 422 F Sequence: 5'-CATCTGCGGACTGGAAAAACAAC-3'
2) Name: 422 Ra Sequence: 5'-GCATCGGTAAACATCTGCTCAAAC-3'

**Primer location**: 422 F and 422 Ra: Reverse tet transactivator gene (rtTA)

Assay Name: MCKrtTA

### **PCR Master Mix Components:**

Component	Manufacturer	Concentration	μl/rxn
Extract-N-Amp PCR Reaction Mix	Sigma (Cat# XNAT2R)	2X	10
422 F	IDT	25µM	0.3
422 Ra	IDT	25µM	0.3
Sterile Water			5.4

#### PCR Setup:

Final Reaction: 16 µl master mix & 4 µl DNA template (10-20 ng/µl)

All reactions were performed in 200µl thin walled PCR tubes and were run in Perkin Elmer 2400 thermocycler or Applied Biosystems 2700 thermocycler.

### **Cycle Parameters:**

1) 94°C 3 minutes 2) 94°C 30 seconds 3) 61°C 30 seconds 4) 72°C 1minute

5) Repeat steps 2-4 34 times for a total of 35 cycles


6) 72°C 10 minutes

7) 4°C hold until refrigerate product

## **Product Analysis:**

All products were analyzed on a 3% agarose gel with ethidium bromide staining Transgene Positive = 402 bp Transgene Negative = no band

## **Example of Gel:**



Lanes 1-9: samples positive for the transgene (402bp band).

Lanes 10 and 11 are extraction and PCR blanks, respectively.

Lane 12 is a WT control (no product).

Lane 13 is a positive control (402bp band).

Lane 14 is 1Kb+ Ladder (Invitrogen Cat# 10787-018)