

09.18.09 MS
07.31.10 HB updated
02.24.14 MLS

Genotyping Protocol: **MMRRC 30406**

Assay Type: PCR- can distinguish heterozygous animals from homozygous animals

DNA Extraction: DNA from tail snips was extracted using Qiagen's DNeasy Blood and Tissue kit (Cat# 69506). Kit directions for animal tissues were performed with a few minor modifications as follows: repeat AW1 and AW2 wash steps one time, elute in 200ul of AE buffer once.

A. PlexinA3 PCR

Strain Description: This strain has IRES-Tau-GFP-Neomycin inserted into exons 2-4 of the Plexin A3 gene (*Plxna3*) on Chromosome X. Details can be found in Yaron et al (2005) Neuron 45:513-523.

Primer Information:

- | | |
|--------------------------|--|
| 1) Name: M30406 141 | Sequence: 5'-CAC CTG CTT CTC ACT CAG GA-3' |
| 2) Name: M30406 143 | Sequence: 5'-GCC AAG TTC TAA TTC CAT CAG A-3' |
| 3) Name: M30406 145 WT B | Sequence: 5'-CCA CAT CCG GCG CCG AAT CCA AT-3' |

Primer location: M30406 141 and 145 WT B are located within the *Plxna3* gene on Chromosome X. M30406 143 is located within the inserted IRES-Tau-GFP-Neomycin.

Assay Name: PlexinA3 PCR

PCR Master Mix Components:

PlexinA3 WT PCR:

component	manufacturer	concentration	μl/rxn
Buffer with MgCl ₂ (green cap)	Roche	10X	2
dNTPs	Promega (Cat# U1515)	1.25mM	3.2
M30406 141	Sigma	25μM	0.3
M30406 145 WT B	Sigma	25μM	0.3
FastStart <i>Taq</i>	Roche (Cat# 12032953001)	5 U/μl	0.2
sterile water			13

PlexinA3 Mut PCR:

component	manufacturer	concentration	μl/rxn
Buffer with MgCl ₂ (green cap)	Roche	10X	2
dNTPs	Promega (Cat# U1515)	1.25mM	3.2
M30406 141	Sigma	25μM	0.3
M30406 143	Sigma	25μM	0.3
FastStart <i>Taq</i>	Roche (Cat# 12032953001)	5 U/μl	0.2
sterile water			13

PCR Setup(both Mut and WT):

Final Reaction: 19μl master mix & 1μl DNA template (10-20ng/μl)

All reactions were performed in 200μl thin walled PCR tubes and were run in Perkin Elmer 2400 thermocycler or Applied Biosystems 2700 thermocycler.

09.18.09 MS
07.31.10 HB updated
02.24.14 MLS

Cycle Parameters (both WT and Mut):

- | | | |
|----|--|--------------------------------|
| 1) | 95°C | 5 minutes |
| 2) | 94°C | 1 minute |
| 3) | 66°C | 1 minute |
| 4) | 72°C | 1 minute |
| 5) | Repeat steps 2-4 34 times for a total of 35 cycles | |
| 6) | 72°C | 10 minutes |
| 7) | 4°C | hold until refrigerate product |

Product Analysis:

Mut PCR products were analyzed on a 3% agarose gel with ethidium bromide staining. WT PCR products were analyzed on a 1% agarose gel with ethidium bromide staining.


Homozygous (female) or Hemizygous (male): 300bp

Heterozygous (female): 1013bp, 300bp

Wild Type: 1013bp


Example gels:

PlexinA3 WT:



Wells 1-10 are negative for the WT allele. Wells 11 and 12 are blanks. Well 13 is a WT control and Well 14 is a homozygous control. Well 15 is 1Kb+ Ladder (Invitrogen Cat# 10787-018).

PlexinA3 Mut:



Wells 1-10 are homozygous female and hemizygous male samples. Wells 11 and 12 are blanks. Well 13 is a WT control and Well 14 is a homozygous control. Well 15: 1Kb+ Ladder.

B. PlexinA4 PCR:

Strain Description: An ER retention signal and Neomycin cassette are inserted into exons 18 and 19 of the Plexin A4 gene (*Plxna4*) on Chromosome 6. Details can be found in Yaron et al (2005) Neuron 45:513-523.

09.18.09 MS
07.31.10 HB updated
02.24.14 MLS

Primer Information:

PlexinA4 PCR:

1) Name: M30406 61	Sequence: 5'-CTT CAG CAC TGG CTG CTG TCA TCT-3'
2) Name: M30406 62	Sequence: 5'-CCA TGC TCT CCT TCA GCC TGC TCT-3'
3) Name: M30406 63	Sequence: 5'-GCT AAA GCG CAT GCT CCA GAC TGC-3'

Primer Location: M30406 61 and 62 are located in the PlexinA4 gene on Chromosome 6. M30406 63 is located in the insertion.

Assay name: PlexinA4 PCR

PCR Master Mix Components:

PlexinA4 PCR:

component	manufacturer	concentration	μl/rxn
Buffer with MgCl ₂ (green cap)	Roche	10X	2
dNTPs	Promega (Cat# U1515)	1.25mM	3.2
M30406 61	Sigma	25μM	0.3
M30406 62	Sigma	25μM	0.3
M30406 63	Sigma	25μM	0.3
FastStart Taq	Roche (Cat# 12032953001)	5 U/μl	0.2
sterile water			12.7

PCR Setup:

Final Reaction: 19μl master mix & 1μl DNA template (10-20ng/μl)

All reactions were performed in 200μl thin walled PCR tubes and were run in Perkin Elmer 2400 thermocycler or Applied Biosystems 2700 thermocycler.

Cycle Parameters:

PlexinA4 PCR:

1)	95°C	5 minutes
2)	94°C	1 minute
3)	68°C	1 minute
4)	72°C	1 minute
5)	Repeat steps 2-4 34 times for a total of 35 cycles	
6)	72°C	10 minutes
7)	4°C	hold until refrigerate product

Product Analysis:

All products were analyzed on a 3% agarose gel with ethidium bromide staining.

PlexinA4 PCR:

Heterozygous: 428bp, 676bp


Homozygous mutant: 428bp

Wild Type: 676bp

09.18.09 MS
07.31.10 HB updated
02.24.14 MLS

Example gel:

PlexinA4:



Wells 1-3, 6-8 are homozygous mutant. Wells 4, 5, 9, 10 are WT. Wells 11 and 12 are blanks. Well 13 is a WT control and Well 14 is a homozygous mutant control. Well 15: 1Kb+ Ladder (Invitrogen Cat# 10787-018).