

11.06.08 MS
08.02.10 HB updated
03.04.14 MLS

Genotyping Protocol: **MMRRC 30410**

Assay Type: PCR- can distinguish heterozygous animals from homozygous animals

DNA Extraction: DNA from tail snips was extracted using Qiagen's DNeasy Blood and Tissue kit (Cat# 69506). Kit directions for animal tissues were performed with a few minor modifications as follows: repeat AW1 and AW2 wash steps one time, elute in 200ul of AE buffer once.

Strain Description: This strain has a secretory gene trap cassette containing beta-gal/neomycin phosphotransferase fusion inserted into the unc-5 homolog B gene (*Unc5h2*) on Chromosome 10. This results in a 2kb deletion which includes exons 3 and 4 of *Unc5h2*. Details can be found in Lu et al (2004) Nature 432:179-186.

Primer Information:

- | | |
|----------------------------|--|
| 1) Name: M30410 Unc5bfwd | Sequence: 5'-ACT AGA ATG CTG TCC AGA C-3' |
| 2) Name: M30410 Unc5brev | Sequence: 5'-AGA GGA GAG CAA CGG ATG-3' |
| 3) Name: M30400 Z2046 anti | Sequence: 5'-TCG TCT GCT CAT CCA TGA CC-3' |
| 4) Name: M30400 Z1924 | Sequence: 5'-GTC TCG TTG CTG CAT AAA CC-3' |

Primer location: M30410 Unc5bfwd and Unc5brev are both located between exons 2 and 3 of the *Unc5h2* gene on Chromosome 10. M30400 Z2046 and Z1924 are located on the inserted gene trap.

Assay name: Unc5h2 PCR

PCR Master Mix Components:

M30410 Mut :

component	manufacturer	concentration	µl/rxn
Buffer with MgCl ₂ (green cap)	Roche	10X	2
dNTPs	Promega (Cat# U1515)	1.25mM	3.2
M30400 Z2046 anti	Sigma or IDT	25µM	0.3
M30400 Z1924	Sigma or IDT	25µM	0.3
FastStart <i>Taq</i>	Roche (Cat# 12032953001)	5 U/µl	0.2
sterile water			13

M30410 WT :

component	manufacturer	concentration	µl/rxn
Buffer with MgCl ₂ (green cap)	Roche	10X	2
dNTPs	Promega (Cat# U1515)	1.25mM	3.2
M30410 Unc5bfwd	Sigma or IDT	25µM	0.3
M30410 Unc5brev	Sigma or IDT	25µM	0.3
FastStart <i>Taq</i>	Roche (Cat# 12032953001)	5 U/µl	0.2
sterile water			13

PCR Setup:

Final Reaction: 19µl master mix & 1µl DNA template (10-20ng/µl)

All reactions were performed in 200µl thin walled PCR tubes and were run in Perkin Elmer 2400 thermocycler or Applied Biosystems 2700 thermocycler.

11.06.08 MS

08.02.10 HB updated

03.04.14 MLS

Cycle Parameters:

Mut PCR:

- 1) 95°C 5 minutes
- 2) 94°C 1 minute
- 3) 68°C 1 minute
- 4) 72°C 1 minute
- 5) Repeat steps 2-4 34 times for a total of 35 cycles
- 6) 72°C 10 minutes
- 7) 4°C hold until refrigerate product

WT PCR:

- 1) 95°C 5 minutes
- 2) 94°C 1 minute
- 3) 63.5°C 1 minute
- 4) 72°C 1 minute
- 5) Repeat steps 2-4 34 times for a total of 35 cycles
- 6) 72°C 10 minutes
- 7) 4°C hold until refrigerate product

Product Analysis:

All products were analyzed on a 3% agarose gel with ethidium bromide staining.

Expected products:


WT: 400bp

Mutant: 500bp


<u>Genotype:</u>	<u>Mut PCR</u>	<u>WT PCR</u>
Heterozygous	500bp	400bp
Homozygous	500bp	no product
Wild Type	no product	400bp

Example gels:

Mut Gel:



WT Gel:



Wells 1 and 2 are positive for the mutant allele (500bp products). Wells 3 and 4 are negative controls. Well 5 is a WT control. Well 6 is a heterozygous control. Well 7 is 1Kb+ Ladder (Invitrogen Cat# 10787-018).

Wells 9 and 10 are positive for the WT allele (400bp products). Wells 11 and 12 are negative controls. Well 13 is a WT control and Well 14 is a heterozygous Well 15 is 1Kb+ Ladder (Invitrogen Cat# 10787-018).

Analysis:

Sample 1 (wells 1 and 9) and Sample 2 (wells 2 and 10) are both Heterozygous.