

FINAL REPORT OF LABORATORY EXAMINATION

4011 Discovery Drive, Columbia, MO 65201

1-800-669-0825 1-573-499-5700

idxxbioanalytics@idexx.com

www.idexxbioanalytics.com

IDEXX BioAnalytics Case # 106474-2025

Received Date: 3/19/2025

Version # 1

Report Date: 3/26/2025

SEROLOGY DETAILS

	Baseline	112-F	112-G	112-H	113-A	113-B	113-C	113-D	113-E	113-F	114-A
LCMV											
LCMV recombinant	MFI > 3.000	-	-	-	-	-	-	-	-	-	-
Mycoplasma pulmonis											
<i>M. pulmonis</i> purified bacteria	MFI > 4.250	-	-	-	-	-	-	-	-	-	-
Pneumocystis carinii											
<i>P. carinii</i>	MFI > 1.500	-	-	-	-	-	-	-	-	-	-
RPV											
RPV purified virus	MFI > 2.500	-	-	-	-	-	-	-	-	-	-
NS1 ¹	MFI > 3.750	-	-	-	-	-	-	-	-	-	-
RMV											
RMV VP2 recombinant	MFI > 2.000	-	-	-	-	-	-	-	-	-	-
NS1 ¹	MFI > 3.750	-	-	-	-	-	-	-	-	-	-
KRV											
KRV purified virus	MFI > 4.250	-	-	-	-	-	-	-	-	-	-
NS1 ¹	MFI > 3.750	-	-	-	-	-	-	-	-	-	-
H1											
H1 purified virus	MFI > 1.750	-	-	-	-	-	-	-	-	-	-
NS1 ¹	MFI > 3.750	-	-	-	-	-	-	-	-	-	-
PVM											
PVM purified virus	MFI > 2.500	-	-	-	-	-	-	-	-	-	-
RCV/SDAV											
RCV/SDAV purified virus	MFI > 3.750	-	-	-	-	-	-	-	-	-	-
RCV/SDAV Spike	MFI > 3.750	-	-	-	-	-	-	-	-	-	-
REO3											
REO3 purified virus	MFI > 4.900	-	-	-	-	-	-	-	-	-	-
RTV											
RTV purified virus	MFI > 2.000	-	-	-	-	-	-	-	-	-	-
TMEV purified virus	MFI > 2.000	-	-	-	-	-	-	-	-	-	-

	Baseline	112-F	112-G	112-H	113-A	113-B	113-C	113-D	113-E	113-F	114-A
Sendai											
Sendai purified virus	MFI > 3.750	-	-	-	-	-	-	-	-	-	-

	Baseline	114-B	114-C	114-D
LCMV				
LCMV recombinant	MFI > 3.000	-	-	-
Mycoplasma pulmonis				
<i>M. pulmonis</i> purified bacteria	MFI > 4.250	-	-	-
Pneumocystis carinii				
<i>P. carinii</i>	MFI > 1.500	-	-	-
RPV				
RPV purified virus	MFI > 2.500	-	-	-
NS1 ¹	MFI > 3.750	-	-	-
RMV				
RMV VP2 recombinant	MFI > 2.000	-	-	-
NS1 ¹	MFI > 3.750	-	-	-
KRV				
KRV purified virus	MFI > 4.250	-	-	-
NS1 ¹	MFI > 3.750	-	-	-
H1				
H1 purified virus	MFI > 1.750	-	-	-
NS1 ¹	MFI > 3.750	-	-	-
PVM				
PVM purified virus	MFI > 2.500	-	-	-
RCV/SDAV				
RCV/SDAV purified virus	MFI > 3.750	-	-	-
RCV/SDAV Spike	MFI > 3.750	-	-	-
REO3				
REO3 purified virus	MFI > 4.900	-	-	-
RTV				
RTV purified virus	MFI > 2.000	-	-	-
TMEV purified virus	MFI > 2.000	-	-	-
Sendai				
Sendai purified virus	MFI > 3.750	-	-	-

NS1¹: NS1 protein is highly conserved among rodent parvoviruses and thus serves as a generic assay for parvovirus seroconversion.

Legend: + = positive - = negative blank = test not performed EQ = equivocal HE = hemolysis precluded testing I = insufficient W = weak positive WB = Western Blot confirmatory analysis pending NS = non-specific reactivity N = normal IgG L = less than normal IgG Positive MFI results are reported as "+" followed by a number from 1 to 33 in thousands rounded off to the nearest thousand.